

Product Nexus Rugged Display

Product Code: NDP-21

Description: Rugged tablet with Nexus software for basic call system. Ideal for rapid deploy nurse call. Must include receiver dongle or ZP-01(M) gateway to operate



Features:

- 10" Industrial-grade Android tablet with docking station for convenient placement on nurse desk
- Receiving options for smaller, medium and large wireless systems
- Smartphone app available with acknowledgement
- Signaling options for remote displays or pagers (requires docking station)

Works with:

- DON-20/30/40 Wireless Receiver Dongle
- ZP-01(M) Master Protektor Unit
- GP2009TR Pager Controller
- NDP-STN Nexus Docking Station

Technical Specifications

Operating frequency: Varies according to dongle
 Operating range: Varies according to dongle
 Power supply: DC5V input, 2000mA
 Battery back-up for 4+ hours
 Dimensions: 280 x 190 x 25mm
 Net weight: 1250g
 Operating temp.: 5°C to 40°C
 Rel. humidity: <90%

Standards

EMC:	EN 301 489-1, 3 &17
	EN 55022,
	EN 55024
	EN 50332-2
	EN 62479
	EN 300 328
Radio:	EN 300 440-1&2
Safety:	EN 60950-1