

Date Order User Address

Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)

<input checked="" type="radio"/> Polished black	22" 24" 25" 26"	
<input type="radio"/> Bright silver	22" 24"	00833
<input type="radio"/> Bright anthracite	24"	00833
<input type="radio"/> Bright orange	22" 24"	00833
<input type="radio"/> Bright fuchsia	22" 24"	00833
<input type="radio"/> Bright blue	22" 24"	00833
<input type="radio"/> Bright red	22" 24"	00833
<input type="radio"/> Bright lime	22" 24"	00833
<input type="radio"/> Pearl white	22" 24"	00833

Rear wheels

<input checked="" type="radio"/> 24" Ultra-light Progeo - Standard!		
<input type="radio"/> 22" Ultra-light Progeo - Standard!		
<input type="radio"/> 25" Ultra-light Progeo - Standard!	(colour: black)	
<input type="radio"/> 26" Progeo - Standard!	(colour: black)	
<input type="radio"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	Weight reduction approx. 450gr.	04383
<input type="radio"/> 24" Ultra-light "Spinergy LX"		01685
<input type="radio"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes <input type="checkbox"/> Hub and spokes black	07846
<input type="radio"/> 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270
<input type="radio"/> 25" Ultra-light "Spinergy LX"		02101
<input type="radio"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes	07847
<input type="radio"/> 26" Ultra-light "Spinergy LX"		02102
<input type="radio"/> 24" Off-road (aluminium handrim black)		01392
<input type="radio"/> 22" with brake handles for companion (drum brake) (aluminium handrim)	Available only with: camber 0° and Setting Standard	03385
<input type="radio"/> 24" with brake handles for companion (drum brake) (aluminium handrim)	Available only with: camber 0° and Setting Standard	03384
<input type="radio"/> One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°	<input type="radio"/> Right <input type="radio"/> Left	02103

Rear wheel tyre

<input checked="" type="radio"/> High pressure tyre	(22"-24"-25"-26")	
<input type="radio"/> High pressure profile tyre	(24")	
<input type="radio"/> Marathon Plus	(24"-25")	with anti-puncture protection 01669
<input type="radio"/> Anti-puncture - Colour Black	(24")	11083
<input type="radio"/> Anti-puncture - Colour Grey	(22")	04181

Handrim **Wheel distance** 3 cm 4 cm


<input checked="" type="radio"/> Aluminium	(22"-24"-25"-26")	
<input type="radio"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 09170
<input type="radio"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 01684
<input type="radio"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 00016
<input type="radio"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 00017
<input type="radio"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 00834
<input type="radio"/> Ergo Para	(24")	00632
<input type="radio"/> TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 11023
Oval countour with a rubber insert		
<input type="radio"/> NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 11025
Ergonomic black profile		
<input type="radio"/> NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers) 11026
Ergonomic profile with rubber coating (ideal for quadriplegics)		
<input type="radio"/> SURGE LT	(25")	Oval countour with a rubber insert 06427
<input type="radio"/> Quad	(24")	00693
<input type="radio"/> Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26") 04050


Wheel locks


<input checked="" type="radio"/> Bent lever push to lock		
<input type="radio"/> Straight lever push to lock		
<input type="radio"/> Straight lever pull to lock		
<input type="radio"/> Sport brake		06428
<input type="radio"/> Brake with extension lever push to lock		01949
<input type="radio"/> Brake with extension lever pull to lock		03372


Date Order User Address


Footplate - See reference tables N°3


- Automatic Closure footplate rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 ,00
  11526

- One piece footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 ,00
  11527

- One piece aluminium footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal 2/3 External External
  11528

- Plastic footrest rake adjustable
 Minimum footplate distance 35cm. 2/3 Internal 2/3 External
 Need a stretch calf band.
NOT AVAILABLE WITH FRONT FORK SUPPORT V DESIGN 2.0
 06861


- One piece Carbon fibre footplate rake adjustable, swing away, flip-up.
 Minimum footplate distance 35cm. Internal External
  03291

- Footplate COMPACT, automatic Closure, rake and depth adjustable
 Minimum footplate distance 35cm. Internal External
  09957


Mudguard


- Mudguard**
- ABS mudguard
 - ABS Straight sideguard
 - Carbon fibre mudguard
 - Carbon fibre Straight sideguard
- 
- 03946
-
- 08918
-
- Removable version**
- Removable ABS mudguard
 - Removable ABS fibre straight sideguard
 - Removable carbon fibre mudguard
 - Removable Carbon fibre straight sideguard
- 
- 08921
-
- 08921
-
- 08919
-
- 08920

Armrests

- Tubular armrests - Steel support with safety locking system
 (Swing away, removable and height adjustable armrest)
Not available with anti-tippers if the rear seat height is less than 40cm
 00694

- Removable armrest with ABS mudguard**

 Desk armrests PROGEO (Removable and height adjustable)
  09956

- Removable armrest with ABS straight sideguard
 Removable armrest with carbon-fiber straight sideguard
  08881
08882





Accessories

- Pad for front frame protectors 07058
- Backrest-seat connection cloth 10241
- Removable storage bag (with quick coupling) 09811



Date Order User Address

Optional 1

Only one optional devices for side can be chosen

<input type="radio"/> Transit wheels	Exclude the installation of : Anti-tippers - Tipping aid - Crutches support		00836
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) <input type="radio"/> Right <input type="radio"/> Left		00023 00024
<input type="radio"/> Tipping aid	<input type="radio"/> Right <input type="radio"/> Left		00025 01046
<input type="radio"/> Crutches support	<input type="radio"/> Right <input type="radio"/> Left		00838 00837

Optional 2

<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947
<input type="radio"/> Quick release front forks	To be matched only with: Standard support front fork		00030
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048
<input type="radio"/> Safety pelvic belt			03972
<input type="radio"/> Harness			04103
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013
<input type="radio"/> Swing away lateral support Not available with Posture backrest PHYSIO	<input type="radio"/> Right <input type="radio"/> Left		02097 02098
<input type="radio"/> Spoke guards Black	(22"-24")		01043
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106

Additional Kit


<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357

Date Order User Address


PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18+26mm
Depth adjustment (4cm)
Angle adjustment ($\pm 18^\circ$)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40			
Shell's shape  7,5 cm	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape  7,5 cm	Lateral retaining wing B=regular		
	R	R	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Configuration DEEP with side support PAD (height adjustable)

- Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



Right
10272



Left
10272

Large (h 20cm)




Right
10273



Left
10273

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



- Large (18x12)

05647

- Small (14x9)

05643

- Lateral support extension

Quantity (number pcs).....

10271

- Lateral support extension

Quantity (number pcs).....

- Small (14x9)

05642

- Large (18x12)

05646

- Headrest adjustable - Double ball mounting system

Available only with backrest Physio
Not available for backrest Physio AIR carbon shell



- Large (34x14)
- Small (28x14)



05648

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo TEKNA ADVANCE - Rev. 2023.2	Prices - VAT excluded to be added to the price list	<p style="text-align: center; margin: 0;">Send this page only if any optional have been chosen</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%; text-align: center;">☆</td> <td style="width:70%;">standard</td> <td style="width:20%; text-align: right; font-size: small;">(supplied if no requests are indicated)</td> </tr> <tr> <td style="text-align: center;">□</td> <td>on request</td> <td style="text-align: right; font-size: small;">(without extra charge)</td> </tr> <tr> <td style="text-align: center;">○</td> <td>optional</td> <td style="text-align: right; font-size: small;">(cost to be added at the basic price)</td> </tr> </table>	☆	standard	(supplied if no requests are indicated)	□	on request	(without extra charge)	○	optional	(cost to be added at the basic price)
☆	standard	(supplied if no requests are indicated)									
□	on request	(without extra charge)									
○	optional	(cost to be added at the basic price)									

REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SMALL 	max 49	max 50	max 51	
MEDIUM 	max 53	max 54	max 55	max 57
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
	Rear wheels			
	22"	24"	25"	26"
H min	min. 36	min. 38	min. 39	min. 40
H max	max 44	max 48	max 49	max 51

3	Footplate option	Footplate distance
	Automatic Closure footplate	min. 36
	Automatic Closure footplate , with multifunction rotative system	min. 31
	One piece footplate - Flip-up!	min. 36
	One piece footplate - Flip-up!, with multifunction rotative system	min. 31
	One piece aluminium footplate - Flip-up!	min. 36
	Plastic footrest	min. 36
	One piece Carbon fibre footplate	min. 36
Calculated with the best configuration and with the footplates positioned on inside.		
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)

Date	Order	User	Address
------	-------	------	---------

12" wheel - Motor 250W (750W maximum power)
14" wheel - Motor 1000W (1500W maximum power)

- Color display with switcher
- Maximum speed limit settable from display
- 5 speed levels selectable
- 10 display-selectable acceleration levels
- Adaptive Cruise Control
- Disc and electronic brake control lever
- Motor brake (adjustable in 10 levels from the display)
- Electronic brake with control bottom (adjustable in 10 levels from the display)
- Power regain when braking
- Error messages management from display (with error history)
- Rear gear beeper
- 48V - 5.8 Ah battery Autonomy up to 25km
- 48V - 11.6 Ah Power plus battery Autonomy up to 50km
- Battery charger 48V
- Rear gear
- USB port
- Easle
- Mudguard
- Front spoiler with predisposition for FULL LED light

(Certified software)

Mototronik

10979




Handlebar

☆ STANDARD controls

☆ Thumb accelerator and Disc and electronic brake control lever (Right)


Grip throttle and Disc and electronic brake control lever



Right
 Left


10913
10914

QUADRIPLÉGIC controls
Handlebar PUSH&PULL




10905

☆ **MOTOTRONIK 12" wheel**
Motor 250W (750W maximum power)
Maximum speed according to law



MOTOTRONIK 14" wheel
Motor 1000W (1500W maximum power)
Maximum speed according to law



10903

Battery

☆ 48V - 5.8 Ah battery *Autonomy up to 25 km*


48V - 11.6 Ah Power plus battery *Autonomy up to 50 km*

10904

Date Order User Address


Accessories

- Handlebar extension
 Adjustable from 15cm up to 23cm



10907

- Motor unit facilitated unlock lever
 Recommended with handlebar for quadriplegic




10915

- Additional FLY 48V - 2.9 Ah battery
 Autonomy up to 15km
 Approved for air transport

10906


- FULL LED light (high intensity)
 Integrated in the spoiler



10908


- Rearview mirror

Right
 Left




10909
10909

- Basket (removable) - (with KLICKfix® attachment system)



10912

- Removable front weight - (with KLICKfix® attachment system)
 (it increases the weight over the power wheel to improve traction)



10910

- Travel bag

10916

Note

Send this page only if any optional have been chosen	
☆	standard (supplied if no requests are indicated)
□	on request (without extra charge)
○	optional (cost to be added at the basic price)