

Date Order User Address

Rear wheels

- ☆ 24" Ultra-light Progeo - Standard!
- 22" Ultra-light Progeo - Standard!
- 25" Ultra-light Progeo - Standard! (colour: black)
- 26" Progeo - Standard! (colour: black)
- 24" ULTRA-LIGHT PROGEO KEY SPORT **Weight reduction approx. 450gr.** 04383
- 24" Ultra-light "Spinergy LX" 01685
- 24" Ultra-light "Spinergy LXL" Black hub and white spokes Hub and spokes black 07846
- 24" Ultra-light "Spinergy Flex rim" **Integrated aluminium/rubber handrim** 02270
- 25" Ultra-light "Spinergy LX" 02101
- 25" Ultra-light "Spinergy LXL" Black hub and white spokes 07847
- 26" Ultra-light "Spinergy LX" 02102
- 24" Off-road (aluminium handrim black) 01392
- 22" with brake handles for companion (drum brake) (aluminium handrim) **Available only with: camber 0° and Setting Standard** 03385

- 24" with brake handles for companion (drum brake) (aluminium handrim) **Available only with: camber 0° and Setting Standard** 03384

- One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°
 - Right Left
 02103





Rear wheel tyre

- ☆ High pressure tyre (22"-24"-25"-26")
- High pressure profile tyre (24")
- Marathon Plus (24"-25") **with anti-puncture protection** 01669
- Anti-puncture - Colour Black (24") 11083
- Anti-puncture - Colour Grey (22") 04181





Handrim **Wheel distance** ☆ 3 cm 4 cm

- ☆ Aluminium (22"-24"-25"-26")
- Aluminium black SOFT-TOUCH (24"-25") (24" handrim, on 25" wheel is installed with spacers) 09170
- Chromed Steel (24"-25") (24" handrim, on 25" wheel is installed with spacers) 01684
- Titanium (24"-25") (24" handrim, on 25" wheel is installed with spacers) 00016
- Ultra Grepp (24"-25") (24" handrim, on 25" wheel is installed with spacers) 00017
- Max Grepp (24"-25") (24" handrim, on 25" wheel is installed with spacers) 00834
- Ergo Para (24") 00632
- TOP FRICTION (24"-25") (24" handrim, on 25" wheel is installed with spacers) 11023
- Oval countour with a rubber insert**
- NOVA H (24"-25") (24" handrim, on 25" wheel is installed with spacers) 11025
- Ergonomic black profile**
- NOVA D (24"-25") (24" handrim, on 25" wheel is installed with spacers) 11026
- Ergonomic profile with rubber coating (ideal for quadriplegics)**
- SURGE LT (25") **Oval countour with a rubber insert** 06427
- Quad (24") 00693
- Silicon cover for standard aluminium handrim Available for measures (22"- 24" - 26") 04050

Wheel locks

- ☆ Bent lever push to lock 
- Straight lever push to lock 
- Straight lever pull to lock 
- Sport brake  06428
- Brake with extension lever push to lock  01949
- Brake with extension lever pull to lock  03372

Footplate - See reference tables N°3

- ☆ Plastic footrest rake adjustable  11525
 - Minimum footplate distance 35cm. ☆ 2/3 Internæ 2/3 External
- One piece footplate - Flip-up! rake and depth adjustable  11527
 - Minimum footplate distance 35cm. Internal ☆ 2/3 Internal External
 - Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 ,00
- One piece aluminium footplate - Flip-up! rake and depth adjustable  11528
 - Minimum footplate distance 35cm. Internal ☆ 2/3 Internal 2/3 External External
- One piece Carbon fibre footplate rake adjustable, swing away, flip-up.  03291
 - Minimum footplate distance 35cm. Internal External

Date Order User Address

Mudguard

Mudguard

- ABS mudguard
- ABS Straight sideguard
- Carbon fibre mudguard
- Carbon fibre Straight sideguard

Removable version

- Removable ABS mudguard
- Removable ABS fibre straight sideguard
- Removable carbon fibre mudguard
- Removable Carbon fibre straight sideguard



03946
08918




08921
08921
08919
08920

Armrests


- Tubular armrests - Steel support with safety locking system
(Swing away, removable and height adjustable armrest)
Not available with anti-tippers if the rear seat height is less than 40cm

Removable armrest with ABS mudguard


- Desk armrests PROGEO (Removable and height adjustable)
- Removable armrest with ABS straight sideguard
- Removable armrest with carbon-fiber straight sideguard




00694



09956



08881



08882

Accessories

- Pad for front frame protectors 07058
- Backrest-seat connection cloth 10241
- Removable storage bag (with quick coupling) 09811

Optional 1


Only one optional devices for side can be chosen

- Transit wheels
- Anti-tipper with rapid hookup
- Tipping aid
- Crutches support


Exclude the installation of :
Anti-tippers -Tipping aid - Crutches support

(inside rotation to reduce the encumbrance)


- Right
- Left
- Right
- Left
- Right
- Left




00836



00023
00024





00025
01046



00838
00837

Optional 2

- Tetraclip (Quick release axle for quadriplegic) 03947
- Quick release front forks To be matched only with: Standard support front fork 00030
- Push handles for companion (Removable and height adjustable) 01048
- Safety pelvic belt 03972
- Harness 04103
- Lumbar support (to insert in the backrest upholstery) 00013
- Swing away lateral support Right 02097
Not available with Posture backrest PHYSIO Left 02098
- Spoke guards Black (22"-24") 01043
- Spoke guards Transparent (22"-24") 00835
- Spoke guards Progeo design (22"-24") 00019
- Spoke guards Dragon (22"-24") 02106

Additional Kit


- 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour) 03378
- 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour) 09357

Date Order User Address


PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18+26mm
Depth adjustment (4cm)
Angle adjustment ($\pm 18^\circ$)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40			
Shell's shape  7.5 cm	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape  7.5 cm	Lateral retaining wing B=regular		
	R	R	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Configuration DEEP with side support PAD (height adjustable)

- Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)



Right
10272

Left
10272

Large (h 20cm)




Right
10273

Left
10273

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



Large (18x12)

05647

Small (14x9)

05643

Lateral support extension

Quantity (number pcs).....

10271

Lateral support extension

Quantity (number pcs).....

Small (14x9)

05642

Large (18x12)

05646

- Headrest adjustable - Double ball mounting system
Available only with backrest Physio
Not available for backrest Physio AIR carbon shell



Large (34x14)
 Small (28x14)




05648

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo TEKNA ADVANCE SWING with footrests removable and Swing away - Rev. 2023.2	Prices - VAT excluded to be added to the price list	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td colspan="2" style="text-align: center; padding: 2px;">Send this page only if any optional have been chosen</td> </tr> <tr> <td style="padding: 2px;"><input checked="" type="checkbox"/> ☆</td> <td style="padding: 2px;">standard (supplied if no requests are indicated)</td> </tr> <tr> <td style="padding: 2px;"><input type="checkbox"/></td> <td style="padding: 2px;">on request (without extra charge)</td> </tr> <tr> <td style="padding: 2px;"><input type="checkbox"/> ○</td> <td style="padding: 2px;">optional (cost to be added at the basic price)</td> </tr> </table>	Send this page only if any optional have been chosen		<input checked="" type="checkbox"/> ☆	standard (supplied if no requests are indicated)	<input type="checkbox"/>	on request (without extra charge)	<input type="checkbox"/> ○	optional (cost to be added at the basic price)
Send this page only if any optional have been chosen										
<input checked="" type="checkbox"/> ☆	standard (supplied if no requests are indicated)									
<input type="checkbox"/>	on request (without extra charge)									
<input type="checkbox"/> ○	optional (cost to be added at the basic price)									

REFERENCE TABLES

1		Front height				
With closed frame (95°) max 5" caster						
FORK	Casters					
	3"	4"	5"	6"		
SMALL 	min. 39 max 42	min. 40 max 43	min. 42 max 43	min. 44 max 45		
MEDIUM 	min. 43 max 45	min. 44 max 46	min. 44 max 47	min. 46 max 48		
LONG 	min. 46 max 52	min. 47 max 53	min. 46 max 54	min. 47 max 55		
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.						

2		Rear height			
		Rear wheels			
		22"	24"	25"	26"
H min		min. 36	min. 38	min. 39	min. 40
H max		max 44	max 48	max 49	max 51

3	Footplate option	Footplate distance
	Elastic footrest.....	min. 36.....
	One piece footplate - Flip-up!, with multifunction rotative system	min. 31
	One piece aluminium footplate - Flip-up!	min. 36
	One piece footplate - Flip-up!	min. 36
	One piece Carbon fibre footplate	min. 36
Calculated with the best configuration and with the footplates positioned on inside.		
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)