



Heating Units

The Relief Pak® moist heat product line just got better!

Meet the newest additions to the Relief Pak® product line: **Relief Pak® Heating Units**

Relief Pak® heating units make using moist heat packs safer and more efficient. Each unit (4, 6, 8, 12, and 24 pack capacity) is made from durable, seamless stainless steel with insulation for both safety and efficiency.

The Relief Pak® heating units feature insulation which truly makes them unique. While the unit is at its set operating temperature, the exterior remains safe to touch. In addition, the insulation helps the heating unit maintain its temperature, which helps to reduce operating costs.

Relief Pak® heating units complement FEI's existing line of heat therapy products. Clinic owners and therapists can trust FEI to provide a complete line of products to suit their every need. Units have a one year warranty. Parts and service are available.



Stationary



4-Pack Unit



6-Pack Unit



8-Pack Unit

Mobile



12-Pack Unit



24-Pack Unit



FABRICATION ENTERPRISES INC
PO Box 1500, White Plains, NY 10602 (USA)
tel: +1-914-345-9300
fax: +1-914-345-9800
FabEnt.com sales@FabEnt.com



AUTHORIZED CE REPRESENTATIVE
AJW Technology Consulting GmbH
Breite Strasse 3
40213 Düsseldorf (Germany)

Relief Pak® Heating Units

- All heating units are insulated (including 4 and 6 pack units)
- Durable, seamless stainless steel for easy cleaning
- Heats moist heat packs in a water bath up to 160°F / 71°C
- On/off switch illuminates green when unit is on
- Heating light glows orange when unit is warming to desired temperature
- Comes with removable hospital-grade plug set



Heating Units Only		
stationary	110 Volt	220 Volt
4 Pack	11-1981	11-1991
6 Pack	11-1982	11-1992
mobile		
8 Pack	11-1984	11-1994
12 Pack	11-1985	11-1995
24 Pack	11-1986	11-1996

Heating Units with Packs

	110 Volt	220 Volt
4 Pack Stationary Units		
4 standard packs	11-1961	11-1971
3 standard, 1 neck pack	11-1961-2	11-1971-2
2 standard, 2 neck packs	11-1961-3	11-1971-3
6 Pack Stationary Units		
3 standard, 1 neck, 2 oversize packs	11-1962	11-1972
6 standard packs	11-1962-2	11-1972-2
8 Pack Mobile Units		
8 standard packs	11-1964	11-1974
4 standard, 4 neck packs	11-1964-2	11-1974-2
2 standard, 2 neck, 2 oversize packs	11-1964-3	11-1974-3
12 Pack Mobile Units		
12 standard packs	11-1965	11-1975
3 standard, 3 neck, 3 oversize packs	11-1965-2	11-1975-2
4 standard, 4 neck, 2 oversize packs	11-1965-3	11-1975-3
4 standard, 4 oversize packs	11-1965-4	11-1975-4
24 Pack Mobile Units		
24 standard packs	11-1966	11-1976
8 oversize, 8 neck packs	11-1966-2	11-1976-2
6 standard, 6 neck, 6 oversize packs	11-1966-3	11-1976-3



Relief Pak® Moist Heat Packs & Covers

- Reusable
- Easy to heat in hot pack heating unit
- Maintains therapeutic moist heat for up to 30 minutes
- Use with Relief Pak® cover or towel
- Washable foam-filled terry cloth or all-terry microfiber covers available for each size



moist heat packs			
	size	each	bulk (12 each)
standard	10" x 12"	11-1310	11-1310-12
neck	7" x 24"	11-1311	11-1311-12
oversize	15" x 24"	11-1312	11-1312-12
half size	5" x 12"	11-1313	11-1313-12
spine small	10" x 18"	11-1314	11-1314-12
spine large	10" x 24"	11-1315	11-1315-12
knee-shoulder	10" x 20"	11-1316	11-1316-12
circular	10" dia	11-1205	11-1205-12

foam-filled hot pack covers			
	size	each	bulk (12 each)
standard	20" x 24"	11-1360	11-1360-12
neck	9" x 24"	11-1361	11-1361-12
oversize	24" x 36"	11-1362	11-1362-12
half size	12" x 15"	11-1363	11-1363-12
circular	12" dia	11-1226	11-1226-12
standard with pocket		11-1364	11-1364-12
oversize with pocket		11-1365	11-1365-12

all-terry covers			
	size	each	bulk (12 each)
standard	20" x 24"	11-1366	11-1366-12
neck	9" x 24"	11-1367	11-1367-12
oversize	24" x 36"	11-1368	11-1368-12
half size	12" x 15"	11-1369	11-1369-12
circular	12" dia	11-1227	11-1227-12

