

Product Pad Transmitter Module

Product Code: TM-11(V3)

Description: Transmitter module for use with Cordless pads and mats. Works with Cordless monitors, Expanders and Pagers.

Features:

- Transmits wireless signals from Cordless bed pads, chair pads, and floor mats to Rondish monitors
- Switch to change between bed & chair type alarms
- Guaranteed for one-year of operation
- Low battery indication

Works with:

- CBM-03 Cordless Safemon
- NGM-21 CordFree Monitor
- ZP-01 Zone Protektor
- CMEX-11/21 Cordless Expander
- CBP-01 Cordless Bed Pad
- CCP-01 Cordless Chair Pad
- CMAT-02 Cordless Floor Mat
- PAG-11/21/31 Caregiver Alert Pager



Technical Specifications

Operating frequency: 433.92MHz
 Modulation type: ASK
 Operating range: 100m (line of sight)
 Power supply: 3V Li ion battery
 Battery life: 1 year
 Dimensions: 7 x 5.5 x 1 cm
 Housing material: ABS
 Net weight: 40g
 Operating temp.: 5°C to 40°C
 Rel. Humidity: < 85%

Standards

EMC: EN 55022
 Radio: EN 300220-1&2, EN 301489-1&3
 Social Alarm: EN 50130-4
 Safety: EN 60950-1
 Medical Device
 Immunity: EN 60601-1-2
 Manufacture: ISO 9001:2008, ISO 13485:2003
 FCC Approved FCCID: WNGTM-11