

Product Universal Under Mattress Bed Pad

Product Code: UMBP-10-U

Description: Sensorpad to trigger an alert when patient exits their bed, works beneath most foam mattresses. USB connection for transmitter or cable

Features:

- Triggers alarm when patient exits their bed
- Usable with either cable or transmitter
- Eyelets for securing to a bed frame. Works best under foam mattresses 14-16cm inches thick
- Durable and long-lasting
- Care instructions written on pad
- 100% latex free

Works with:

- TM-11/12U Cordless Transmitter
- USBC-2M 92" cable, USB to RJ10
- RBC-02 Rubber Band Clip Set
- NGM-11 Palomino Monitor (wired)
- NGM-21 CordFree Monitor (wireless)
- CMEX-11/21 Cordless Expander
- PAG-11/21/31 Caregiver Alert Pager
- ZP-01 Zone Protektor



Technical Specifications

Cable: 2m with RJ10 plug (NO) and USB break point

Material: PVC & PET

Dimensions: 83 x 26cm

Net weight: 620g

Operating temp.: 5°C to 40°C

Operating weight: 30 - 200kg

Internal resistance: < 1200 Ω

Recommended cleaning: 1:100 solution of household bleach (5.25 - 6.15% sodium hypochlorite)

Standards

EMC

Emmissions: EN 55022:2010

Social Alarms: EN 50130-4:1995

Manufacture: ISO 9001:2015, ISO 13485:2016